

Disposable Laryngeal Mask Connector



TRI-1

- Lower-End port delimited for SafeLM® Video Laryngeal Mask System
- Standard 15 mm Horizontal Port for anesthetic breathing circuit.
- Capped Luer-Locked Port on top provides protective access
- Suitable for various catheters such as secretion suction tube.



TRI-2

- One standard 15 mm male connector and one 15 mm female connector for connecting Video Laryngeal Mask with anesthetic breathing circuit.
- Silicone double flip cap with Access Port
- Double flip cap allow passage of bronchoscope. Access Port allow passage for secretion suction tube
- Allows flexible position adjustment of laryngeal mask and breathing circuit while minimizing torque.

Used with Video Laryngeal Mask System



Videoscope



Disposable Video Laryngeal Mask